

#654

PULSE OXIMETER FINGER WRAP



Benefits:

- Pulse oximeter probe wrap that helps secure various styles of optical transducers and probe sensors
- Prevents the misalignment of optical components while blocking out extra light
- Fits a neonatal foot, pediatric toes, or pediatric and adult fingers
- Packaged 12 wraps per package
- Velcro® tab eliminates need for adhesive tape
- Skin friendly foam laminate made from moisture repellent material
- Over-the-counter intended for non-medical, recreational use
- Adjustable, one size fits all
- Secures Pulse Ox sensor to allow for steady reading
- Measures 5" long x 1" wide and can be cut to fit any length