


TRACH-TIE® Tracheostomy Tube Neckband

Product #501

- 100% Cotton Flannel and Open Loop with "Moisture Wicking" Fabric
- Two Piece Adjustable Neckband for Universal Sizing
- Thicker Foam Padding



- "No Roping" Effect with 100% Cotton Non-Stretch Fabric
- Closed Seams on Edge of Flannel Neckband
- "Flex-Piece" for Cough Reflex
- Latex Free

PEPPER MEDICAL

Providing Superior Products at Substantial Savings!

P.O. Box 30914 • Bethesda, MD 20824

1-800-647-0172

www.peppermedical.com



TRACH-TIE® Tracheostomy Tube Neckband

Product #501



CLINICAL ADVANTAGES

Features	Benefits	Clinical Advantages Compared to Others
100% cotton flannel and open loop with "Moisture Wicking" properties-	<ul style="list-style-type: none"> • Soft, Hypo-allergenic • Wicks moisture away from skin through cotton flannel to be evaporated through loop backing 	<ul style="list-style-type: none"> • Prevents irritation and skin breakdown • Keeps skin drier and cooler • Air easily moves through to allow skin to breath
"No Roping" Cotton Fabric	<ul style="list-style-type: none"> • Fabric integrity will not stretch, pull or "rope" digging into neck • 1" width will not roll-over or fold 	<ul style="list-style-type: none"> • "Roping effect" causes skin irritation and pain. • Prevents digging into neck
Latex – Free	<ul style="list-style-type: none"> • No Latex Reactions 	<ul style="list-style-type: none"> • Prevents Latex Hypersensitivity
Two Piece Adjustable neckband for universal sizing	<ul style="list-style-type: none"> • Easy and Quick to use • Easily Adjusted to any size patient • Secure Velcro® Tabs 	<ul style="list-style-type: none"> • Reduces Staff time • Increase patient care and ease of application
"Closed Seams" on edge of flannel Neckband	<ul style="list-style-type: none"> • Reduces sponge effect of foam laminate • Keeps moisture out of trach tie 	<ul style="list-style-type: none"> • Prevents moisture retention & pooling • Reduces skin irritation
Thicker Foam Padding	<ul style="list-style-type: none"> • More Comfortable Patient 	<ul style="list-style-type: none"> • More Comfort- Less Pain
"FlexPiece" for cough reflex	<ul style="list-style-type: none"> • Allows for edema & cough reflex • Allows "Flex" for weight of Vent Circuitry 	<ul style="list-style-type: none"> • Adds "Flex" for patient comfort and safety
Packaged in Boxes of 20	<ul style="list-style-type: none"> • Less boxes stocked • Less Material Management Time 	<ul style="list-style-type: none"> • LESS EXPENSIVE