



210 Carter Drive Unit 10, West Chester, PA 19382
Phone 800-647-0172 Fax 610-692-8468 www.peppermedical.com

Pleated HEPA Filter

#HEPA701



Product Code	HEPA701
Description	Pleated HEPA Filter
Humidity Output	25mg H ₂ O @ VT 500ml
Resistance @ 30 LPM	150pa or 1.53cm H ₂ O
Dead space	52ml
Tidal Volume Range	200 – 1500 ml
Filtration Efficiency	BFE 99.9999% VFE 99.9999%
Weight	40g
Dimensions	h-81.5 mm w- 68.5 mm
Materials of Construction	Filter Media: Hydrophobic Media Housing: Polypropylene
Gas Monitoring Luer Port	Yes
Connections	22M-15F / 22F-15M

The Pepper Medical HEPA Filter come with efficiencies of 99.999%, with luer ports, and an integral cap-and-strap as an additional safety feature to remove the risk of a lost cap. Pleated media filters provide higher efficiencies and greater levels of hydrophobicity.