



600 W Strasburg Road, West Chester, PA 19382  
Phone 800-647-0172 Fax 610-692-8468 www.peppermedical.com

## Heat and Moisture Exchanger #HME001



Product Code	HME001
Description	Tracheal HME
Humidity Output	23.4mg @ VT 500ml
Resistance @ 30 LPM	40pa
Dead space	11ml
Tidal Volume Range	> 25ml

The Pepper Medical Tracheal Heat & Moisture Exchanger captures heat and moisture and returns it to the patient. It is meant for use by spontaneously breathing patients with a tracheostomy tube.

### Product Features

- Low profile and Lightweight.
- High performance paper element.
- Warms and humidifies air with minimal resistance to flow.
- Convenient oxygen delivery.
- Full compatibility with breathing systems.

### Instructions

- To be used on spontaneously breathing patients with a trach tube
- Should be connected to the 15/22 mm connection of the trach tube
- Device should be used for maximum of 24 hours