



Closed Suction Catheter



Features and Benefits:

- **Safe and High-Efficiency:** Secure swivel & T-style fittings offer protection for patients and medical staff, reducing the risk of nosocomial infection and proper ventilation. The suction catheter tip is soft & smooth reducing damage to the airway mucosa after repeated suctioning.
- **Easy Operation Features include:** “Easy-Switch” for Ergonomic suction switch design, with locking device; clear and transparent protective cover for easy observation; double swivel connector design to reduce torque for improved patient comfort. Equipped with extendable flex adapter and separators for easy connection and separation.
- **Full Line:** For adults and children; 24 hours and 72 hours.
- Conforms to the international color code standards to distinguish sizing; **equipped with a color coded Day sticker.**

Specifications:

| PEPPER MEDICAL ITEM# | DESCRIPTION |
|----------------------|---|
| CS24-6 | 6Fr Elbow Closed Suction Catheter |
| CS24-8 | 8Fr Elbow Closed Suction Catheter |
| CS24-10 | 10Fr Elbow Closed Suction Catheter |
| CS24-DS-12 | 12Fr Elbow with Double Swivel Closed Suction Catheter |
| CS24-DS-14 | 14Fr Elbow with Double Swivel Closed Suction Catheter |
| CS24-T-14 | 14Fr T-Piece Closed Suction Catheter |
| CS24-DS-16 | 14Fr Elbow with Double Swivel Closed Suction Catheter |

All of the above available in Trach Tube and Endo Trach Tube lengths.

ORDERING INFORMATION:

1-800-647-0172 | 610-692-8468 | Faxorders@peppermedical.com | www.PepperMedical.com