

CLOSED SUCTION CATHETERS



Features and Benefits:

- **Safe and High-Efficiency:** Secure swivel & T-style fittings offer protection for patients and medical staff, reducing the risk of nosocomial infection and proper ventilation. The suction catheter tip is soft & smooth reducing damage to the airway mucosa after repeated suctioning.
- **Easy Operation Features include:** “Easy-Switch” for Ergonomic suction switch design, with locking device; clear and transparent protective cover for easy observation; double swivel connector design to reduce torque for improved patient comfort. Equipped with extendable flex adapter and separators for easy connection and separation.
- **Full Line:** For adults and children; 24 hours and 72 hours.
- Conforms to the international color code standards to distinguish sizing; **equipped with a color coded Day sticker.**

Configuration:


Conventional configuration:

- 1 Extendable tube
- 2 Separator
- 3 Color code card (Date)
- 4 Endotracheal tube conversion (for Children)
- 5 Medicine injection kit (for Children)



CLOSED SUCTION CATHETERS

Specifications:

 <p>H1 (Normal Type) 24 Hours</p> <p>H2 (Valve Type) 72 Hours</p>	PRODUCT CODE	DESCRIPTION	LENGTH
	CS24-6	Closed Suction Catheter, with Elbow Adaptor, 6Fr	30cm
	CS24-8	Closed Suction Catheter, with Elbow Adaptor, 8Fr	20cm
	CS24-10	Closed Suction Catheter, with Elbow Adaptor, 10Fr	20cm
	CS24-DS-12	Closed Suction Catheter, with Double Swivel Elbow, 12Fr	30cm
	CS24-DS-14	Closed Suction Catheter, with Double Swivel Elbow, 14Fr	30cm
	CS24-T-14	Closed Suction Catheter, with T-Piece, 14Fr	30cm
	CS24-T-16	Closed Suction Catheter, with T-Piece, 16Fr	60cm

**All of the above available in both 24 and 72 hour versions*

ORDERING INFORMATION:

1-800-647-0172 | 610-692-8468 Fax

orders@peppermedical.com | www.PepperMedical.com