



PRODUCT CATALOG

**PEPPER MEDICAL**

1-800-647-0172 | [PepperMedical.com](http://PepperMedical.com)

# TABLE OF CONTENTS



This is an interactive pdf. Hover and click for links to other pages and website.

## CLOSED SUCTION CATHETERS

Closed Suction Catheters **4**

## ENDO-TIES® AND HOLDERS

#790 **5**  
#790B **6**  
#777 **7**

## TRACH-TIES® - Tracheostomy Tube Neckbands

PRODUCT LIST **8**  
#300 **9**  
#301, #301P **10 - 11**  
#501, #501P & #501W **12 - 14**  
#502, #502P **15 - 16**  
#503, #504 **17**  
#540, #541 & #542 **18**  
#597S, #597M & #597L **19**  
#1001 **20**  
#301PRT, #501PRT **21**

## VENT-TIES® - Anti-Disconnect Devices\*

#401, #401P **22 - 23**  
#402 **24**  
#403 **25**  
#405, #405W **26 - 27**  
#406B **28**  
#497S, #497M & #497L **29**

## FOLEY-TIES® - Foley Catheter Holders

#601 **30**  
#606, #606XL **31 - 32**  
#5445 **33**  
#650 **34**  
#5447 **35**

## HME's - Heat Moisture Exchangers & Filters

PRODUCT LIST **36**  
HME-01P, HME-02P **37 - 38**  
HME-02-Sux **39**  
HME-01F, HME-02F **40 - 41**  
HME-03P, HME-04P **42 - 43**  
HME-Mini-P **44**  
HEPA701 **45**  
HME-711-Mini **46**  
HME-712, HME-712-Mini **47 - 48**  
HME-712-Flex **49**  
BV001, BV002 **50 - 51**

## DISPOSABLE RESPIRATORY DEVICES

Respiratory Masks **52**  
Oxygen Tubing **53**  
Breathing Circuit Tubing **54**  
Pulse Oximeter Finger Wrap **55**

(CONTINUED ON THE NEXT PAGE)

# TABLE OF CONTENTS

## TRACHEOSTOMY SPEAKING VALVES

#710, #711 & #740 **56**

## UNI-BLADE LARYNGOSCOPE

Uni-Blade **57**

## VCS/NG SECUREMENT DEVICES

#701VCS, #701VCS-Pediatric **58**

#701VCS/NG **59**

## AB BINDERS AND AB BELTS

AB Binders **60**

AB Belts **61**

### **\*FDA, CDC, ECRI and JOINT COMMISSION – Reference the dangers of Ventilator Disconnections:**

- FDA- Safety: The FDA reports that approximately 200 deaths and numerous injuries occur annually due to accidental ventilator disconnections.
- CDC- VAP: CDC states unintentional ventilator disconnections are a cause of VAP – Ventilator Acquired Pneumonia.
- JCAHO-& AHA Commercially Available Device: Joint Commission and AHA recommends that facilities utilize a commercially available device, not improvised devices such as rubberbands or tape.
- ALARMS: ECRI rated Alarm Management as the #1 most hazardous health technology issue in 2012. "We routinely see deaths and injuries associated with 'alarm fatigue' Alarms are unsettling to family members, distracting to staff and become dangerous 'white noise' that can mask a real event!!"



## Closed Suction Catheter



### Features and Benefits:

- **Safe and High-Efficiency:** Secure swivel & T-style fittings offer protection for patients and medical staff, reducing the risk of nosocomial infection and proper ventilation. The suction catheter tip is soft & smooth reducing damage to the airway mucosa after repeated suctioning.
- **Easy Operation Features include:** “Easy-Switch” for Ergonomic suction switch design, with locking device; clear and transparent protective cover for easy observation; double swivel connector design to reduce torque for improved patient comfort. Equipped with extendable flex adapter and separators for easy connection and separation.
- **Full Line:** For adults and children; 24 hours and 72 hours.
- Conforms to the international color code standards to distinguish sizing; **equipped with a color coded Day sticker.**

### Specifications:

PEPPER MEDICAL ITEM#	DESCRIPTION
CS24-6	6Fr Elbow Closed Suction Catheter
CS24-8	8Fr Elbow Closed Suction Catheter
CS24-10	10Fr Elbow Closed Suction Catheter
CS24-DS-12	12Fr Elbow with Double Swivel Closed Suction Catheter
CS24-DS-14	14Fr Elbow with Double Swivel Closed Suction Catheter
CS24-T-14	14Fr T-Piece Closed Suction Catheter
CS24-DS-16	14Fr Elbow with Double Swivel Closed Suction Catheter

**All of the above available in Trach Tube and Endo Trach Tube lengths.**

#### ORDERING INFORMATION:

1-800-647-0172 | 610-692-8468 | Faxorders@peppermedical.com | www.PepperMedical.com

#790

**PREVENT®**  
**Endo Trach Tube Holder**



Features & Benefits:

- Slide Rail allows for movement and rotation of the tube to prevent sores and skin breakdown
- Larger strap release Tab allows for easier and quicker tube disengagement
- Medical grade adhesive tape
- Cushioned foam pad relieves pressure and will NOT absorb water or adhere to hair
- Soft foam cushion padding protects the teeth, lips and skin
- Soft adjustable neckband
- Accommodates multiple-size ET tubes
- Less packaging waste

**ORDERING INFORMATION:**

1-800-647-0172 | 610-692-8468 Fax

orders@peppermedical.com | www.PepperMedical.com

**#790B**

**PREVENT®**

**Endo Trach Tube Holder  
with attached Bite Block**



## Features & Benefits:

- Slide Rail allows for movement and rotation of the tube to prevent sores and skin breakdown
- Larger strap release Tab allows for easier and quicker tube disengagement
- Medical grade adhesive tape
- Cushioned foam pad relieves pressure and will NOT absorb water or adhere to hair
- Soft foam cushion padding protects the teeth, lips and skin
- Soft adjustable neckband
- Accommodates multiple-size ET tubes
- Less packaging waste

### ORDERING INFORMATION:

1-800-647-0172 | 610-692-8468 Fax

orders@peppermedical.com | www.PepperMedical.com

#777

## ENDO TUBE HOLDER



### Benefits:

- Secures Endotracheal Tubes, LMA Tubes, and other supraglottic airway tubes.
- Soft foam padding with hook and loop that secures around the back of the neck.
- Opening that allows for suctioning
- Leaves the tube visible for inspection and monitoring
- Screw fastener that allows for easy securement of the tube
- Disposable, for single use

### ORDERING INFORMATION:

1-800-647-0172 | 610-692-8468 Fax

orders@peppermedical.com | www.PepperMedical.com

# TRACH-TIES®

## Tracheostomy Tube

### Neckbands

ITEM#	DESCRIPTION
300	One Piece Adult Trach Tie
301	One Piece Adult Trach Tie, cotton
301P	One Piece Pedi Trach Tie, cotton
501	Two Piece Adult Trach Tie, cotton
501P	Two Piece Pedi Trach Tie, cotton
501W	Bariatric Trach Tie, moisture wicking
502	Two Piece Adult Trach Tie with taper ends, cotton
502P	Two Piece Pedi Trach Tie with taper ends, cotton
503	Two Piece Pedi Trach Tie, cotton
504	Two Piece Adult Trach Tie, cotton
540	Two Piece Adult Trach Tie
541	Two Piece Pedi Trach Tie
542	Two Piece Pedi/Neo Trach Tie
597S	Small Trach Tie
597M	Medium Trach Tie
597L	Large Trach Tie
1001	Two Piece Adult Trach Tie
301PRV	Rainbow Trach Tie, 1 Piece, moisture wicking
501PRV	Rainbow Trach Tie, 2 Piece, moisture wicking

#### ORDERING INFORMATION:

1-800-647-0172 | 610-692-8468 Fax

orders@peppermedical.com | www.PepperMedical.com

# #300

## TRACH-TIE® Adult Tracheostomy Tube Neckband - 1 Piece



### Features & Benefits:

- **LATEX FREE**
- Soft, moisture repellent, fabric ensures patient safety and comfort
- Soft, plush backing on hook tabs provides additional comfort
- Reduces skin breakdown
- Easy to adjust for proper fit up to a 21" neck
- "Closed Seams" on edge of neckband reduces sponge effect of foam laminate and keeps moisture out of trach tie
- Thicker foam padding
- Reimbursement Code A7526

### DIMENSIONS

Length: 17" / Width: 1"

### ORDERING INFORMATION:

1-800-647-0172 | 610-692-8468 Fax

orders@peppermedical.com | www.PepperMedical.com

# #301

## TRACH-TIE® Adult Tracheostomy Tube Neckband - 1 Piece



### Features & Benefits:

- **LATEX FREE**
- Soft, 100% hypoallergenic, cotton fabric ensures patient safety and comfort
- Reduces skin breakdown
- Easy to adjust for proper fit up to a 21" neck
- "Closed Seams" on edge of neckband reduces sponge effect of foam laminate and keeps moisture out of trach tie
- Thicker foam padding
- Reimbursement Code A7526

### DIMENSIONS

Length: 17" / Width: 1"

### ORDERING INFORMATION:

1-800-647-0172 | 610-692-8468 Fax

orders@peppermedical.com | www.PepperMedical.com

# #301P

## PEDI-TIE®

### Pediatric Tracheostomy Tube Neckband - 1 Piece



### Features & Benefits:

- **LATEX FREE**
- Soft, 100% hypoallergenic, cotton fabric ensures patient safety and comfort
- Reduces skin breakdown
- **#301P is made specifically for the Pediatric Patient**
- “Closed Seams” on edge of neckband reduces sponge effect of foam laminate and keeps moisture out of trach tie
- Thicker foam padding
- Reimbursement Code A7526

### DIMENSIONS

Length: 11.5” / Width: .75”

### ORDERING INFORMATION:

1-800-647-0172 | 610-692-8468 Fax

orders@peppermedical.com | www.PepperMedical.com

# #501

## TRACH-TIE® Adult Tracheostomy Tube Neckband - 2 Piece



### Features & Benefits:

- **LATEX FREE**
- Soft, 100% hypoallergenic, cotton fabric ensures patient safety and comfort
- Reduces skin breakdown
- Two piece adjustable neckband for universal sizing
- “Closed Seams” on edge of neckband reduces sponge effect of foam laminate and keeps moisture out of trach tie
- “Flex-Piece” for cough reflex
- Thicker foam padding
- Reimbursement Code A7526

### DIMENSIONS

Length: 22.5” / Width: 1”

### ORDERING INFORMATION:

1-800-647-0172 | 610-692-8468 Fax  
orders@peppermedical.com | www.PepperMedical.com

# #501P

## **PEDI-TIE II®** **Pediatric Tracheostomy Tube** **Neckband - 2 Piece**



### Features & Benefits:

- **LATEX FREE**
- Soft, 100% hypoallergenic, cotton fabric ensures patient safety and comfort
- Reduces skin breakdown
- **#501P is made specifically for the Pediatric Patient**
- Two piece adjustable neckband for universal sizing
- “Closed Seams” on edge of neckband reduces sponge effect of foam laminate and keeps moisture out of trach tie
- Thicker foam padding
- Reimbursement Code A7526

### **DIMENSIONS**

Length: 14” / Width: .75”

### **ORDERING INFORMATION:**

1-800-647-0172 | 610-692-8468 Fax  
orders@peppermedical.com | www.PepperMedical.com

# #501W

## TRACH-TIE® Bariatric Tracheostomy Tube Neckband - 2 Piece



### Features & Benefits:

- **LATEX FREE**
- 1.5" WIDE design prevents roping effect
- Two piece adjustable neckband allows for easy adjustment
- Soft, plush backing on hook tabs provides additional comfort
- Made of DriWick®, a proprietary material that pulls moisture away for the patients skin, reducing skin breakdown and irritation
- Extra long for bariatric patients
- Reimbursement Code A7526

### DIMENSIONS

Length: 28" / Width: 1.5"

### ORDERING INFORMATION:

1-800-647-0172 | 610-692-8468 Fax

orders@peppermedical.com | www.PepperMedical.com

# #502

## TRACH-TIE® Adult Tracheostomy Tube Neckband - 2 Piece



### Features & Benefits:

- **LATEX FREE**
- Soft, 100% hypoallergenic, cotton fabric ensures patient safety and comfort
- 2.25" Velcro® hook tabs, tapered to prevent skin breakdown
- Two piece adjustable neckband for universal sizing
- 0.25" thick foam with moisture repellent backing
- Reimbursement Code A7526

### DIMENSIONS

Length: 18.5" / Width: 1"

### ORDERING INFORMATION:

1-800-647-0172 | 610-692-8468 Fax

orders@peppermedical.com | www.PepperMedical.com

# #502P

## PEDI-TIE®

### Pediatric Tracheostomy Tube Neckband - 2 Piece



## Features & Benefits:

- **LATEX FREE**
- Soft, 100% hypoallergenic, cotton fabric ensures patient safety and comfort
- 2.25" Velcro® hook tabs, tapered to prevent skin breakdown
- Two piece adjustable neckband for universal sizing
- 0.25" thick foam with moisture repellent backing
- Reimbursement Code A7526

### DIMENSIONS

Length: 11" / Width: .75"


### ORDERING INFORMATION:

1-800-647-0172 | 610-692-8468 Fax

orders@peppermedical.com | www.PepperMedical.com

# #503

## TRACH-TIE® Pediatric Tracheostomy Tube Neckband - 2 Piece



### Features & Benefits:

- **LATEX FREE**
- 0.25" thick white foam with 100% cotton laminate
- Reduces skin breakdown
- **#503 is made specifically for the Pediatric Patient**
- Two piece adjustable neckband for universal sizing
- 2 .875" Velcro® hook tabs, stitched
- Reimbursement Code A7526

**DIMENSIONS:** Length: 12" / Width: 1"

# #504

## TRACH-TIE® Adult Tracheostomy Tube Neckband - 2 Piece



### Features & Benefits:

- **LATEX FREE**
- 0.25" thick white foam with 100% cotton laminate
- Reduces skin breakdown
- Two piece adjustable neckband for universal sizing
- 2 .875" Velcro® hook tabs, stitched
- Reimbursement Code A7526

**DIMENSIONS:** Length: 18.5" / Width: 1"

**ORDERING INFORMATION:** 1-800-647-0172 | 610-692-8468 Fax  
orders@peppermedical.com | www.PepperMedical.com


## #540

### TRACH-TIE® Adult Tracheostomy Tube Neckband - 2 Piece

#### Features & Benefits:

- **LATEX FREE**
- 0.25" thick blue foam with moisture repellent backing
- 2 .875" Velcro® hook tabs, stitched
- Reimbursement Code A7526

**DIMENSIONS:** Length: 18.5" (including 2.5" of elastic) / Width: 1"




## #541

### TRACH-TIE® PEDIATRIC Tracheostomy Tube Neckband - 2 Piece

#### Features & Benefits:

- **LATEX FREE**
- 0.1875" thick blue foam with moisture repellent backing
- 2 .875" Velcro® hook tabs, welded
- Reimbursement Code A7526

**DIMENSIONS:** Length: 17.5" / Width: 1"




## #542

### TRACH-TIE® PEDIATRIC / NEONATAL Tracheostomy Tube Neckband - 2 Piece

#### Features & Benefits:

- **LATEX FREE**
- 0.1875" thick blue foam with moisture repellent backing
- 2 .875" Velcro® hook tabs, welded
- Reimbursement Code A7526

**DIMENSIONS:** Length: 10" / Width: 0.75"



**ORDERING INFORMATION:** 1-800-647-0172 | 610-692-8468 Fax  
orders@peppermedical.com | www.PepperMedical.com


## #597S

### TRACH-TIE® Small Tracheostomy Tube Neckband - 1 Piece

#### Features & Benefits:

- **LATEX FREE**
- 0.25" thick blue foam with moisture repellent backing
- 2 .875" x 0.1875 Velcro® hook tabs
- Reimbursement Code A7526

**DIMENSIONS:** Length: 9.5" / Width: 1.25"




## #597M

### TRACH-TIE® Medium Tracheostomy Tube Neckband - 1 Piece

#### Features & Benefits:

- **LATEX FREE**
- 0.25" thick blue foam with moisture repellent backing
- 2 .875" x 0.25 Velcro® hook tabs
- Reimbursement Code A7526

**DIMENSIONS:** Length: 18.5" / Width: 1.25"




## #597L

### TRACH-TIE® Large Tracheostomy Tube Neckband - 1 Piece

#### Features & Benefits:

- **LATEX FREE**
- 0.25" thick blue foam with moisture repellent backing
- 2 .875" x 0.25 Velcro® hook tabs
- Reimbursement Code A7526


**DIMENSIONS:** Length: 23.5" / Width: 1.25"



**ORDERING INFORMATION:** 1-800-647-0172 | 610-692-8468 Fax  
orders@peppermedical.com | www.PepperMedical.com

# #1001

## TRACH-TIE® Adult Tracheostomy Tube Neckband - 2 Piece



### Features & Benefits:

- **LATEX FREE**
- 2.875" Velcro® hook tabs, stitched
- Two piece adjustable neckband for universal sizing
- 0.25" thick blue foam with moisture repellent backing
- Reimbursement Code A7526

### DIMENSIONS

Length: 18.5" (including 2.5" of elastic) / Width: 1"

### ORDERING INFORMATION:

1-800-647-0172 | 610-692-8468 Fax

orders@peppermedical.com | www.PepperMedical.com

# #301PRT

## RAINBOW TIES® Pediatric/Neonatal Tracheostomy Tube Neckband - 1 Piece



### Features & Benefits:

- 100% hypoallergenic, cotton fabric ensures patient safety and comfort
- **Packaged 10 per bag with 3 different color schemes: Variety, Pastel, & Bold**
- **LATEX FREE**
- Made of Orthowick® which pulls moisture away from skin aiding in the reduction of skin breakdown and infection
- Soft Velcro® tabs and plush backing
- Reimbursement Code A7526

**DIMENSIONS:** Length: 11.5" / Width: 0.75"

# #501PRT

## RAINBOW TIES® Pediatric/Neonatal Tracheostomy Tube Neckband - 2 Piece



### Features & Benefits:

- 100% hypoallergenic, cotton fabric ensures patient safety and comfort
- **Packaged 10 per bag with 3 different color schemes: Variety, Pastel, & Bold**
- **LATEX FREE**
- Made of Orthowick® which pulls moisture away from skin aiding in the reduction of skin breakdown and infection
- Soft Velcro® tabs and plush backing
- Reimbursement Code A7526

**DIMENSIONS:** Length: 14" / Width: 0.75"

**ORDERING INFORMATION:** 1-800-647-0172 | 610-692-8468 Fax  
orders@peppermedical.com | www.PepperMedical.com

# #401

## **VENT-TIE® Adult Trach Tube Holder with Ventilator Anti-Disconnect Device - 1 Piece**



### Features & Benefits:

- Provides trach tube security while also securing the ventilator circuitry
- Lifts and supports at the 15 mm site, removing a significant amount of pressure on the distal end of the trach tube – intended to reduce back tracheal wall irritation and initiate cough reflex
- Stabilizes the interface – intended to reduce patient discomfort, movement, and subsequent pulling of the trach, potentially reducing decannulation
- Replaces rubberbands, shoe strings, and other ad hoc securement devices
- **LATEX FREE** - No latex reactions which prevents latex hypersensitivity
- Compatible with close suction, T-pieces elbow connectors, OptiFlow and all trach tubes
- 100% hypoallergenic cotton
- Quick release in case of emergency
- Reimbursement Code A7526

### **DIMENSIONS**

Length: 17" / Width: 1"

### **ORDERING INFORMATION:**

1-800-647-0172 | 610-692-8468 Fax

orders@peppermedical.com | www.PepperMedical.com

# #401P

## **PEDI-VENT-TIE® Pediatric Trach Tube Holder with Ventilator Anti-Disconnect Device - 1 Piece**



### Features & Benefits:

- Provides trach tube security while also securing the ventilator circuitry
- Lifts and supports at the 15 mm site, removing a significant amount of pressure on the distal end of the trach tube – intended to reduce back tracheal wall irritation and initiate cough reflex
- Stabilizes the interface – intended to reduce patient discomfort, movement, and subsequent pulling of the trach, potentially reducing decannulation
- Replaces rubberbands, shoe strings, and other ad hoc securement devices
- **LATEX FREE** - No latex reactions which prevents latex hypersensitivity
- Compatible with close suction, T-pieces elbow connectors, OptiFlow and all trach tubes
- **#401P is made specifically for the Pediatric Patient**
- 100% hypoallergenic cotton
- Quick release in case of emergency
- Reimbursement Code A7526

### **DIMENSIONS**

Length: 11.5" / Width: .75"

### **ORDERING INFORMATION:**

1-800-647-0172 | 610-692-8468 Fax

orders@peppermedical.com | www.PepperMedical.com

# #402

## **VENT-TIE® Adult Trach Tube Holder with Ventilator Anti-Disconnect Device - 2 Piece**



### Features & Benefits:

- Provides trach tube security while also securing the ventilator circuitry
- Lifts and supports at the 15 mm site, removing a significant amount of pressure on the distal end of the trach tube – intended to reduce back tracheal wall irritation and initiate cough reflex
- Stabilizes the interface – intended to reduce patient discomfort, movement, and subsequent pulling of the trach, potentially reducing decannulation
- Replaces rubberbands, shoe strings, and other ad hoc securement devices
- **LATEX FREE** - No latex reactions which prevents latex hypersensitivity
- Compatible with close suction, T-pieces elbow connectors, OptiFlow and all trach tubes
- 100% hypoallergenic cotton with elastic stretch piece
- Quick release in case of emergency
- Reimbursement Code A7526

### **DIMENSIONS**

Length: 22" / Width: 1"

### **ORDERING INFORMATION:**

1-800-647-0172 | 610-692-8468 Fax

orders@peppermedical.com | www.PepperMedical.com

# #403

## **VENT-TIE® Adult Trach Tube Holder with Ventilator Anti-Disconnect Straps (for T-Piece) - 2 Piece**



### Features & Benefits:

- Provides trach tube security while also securing the ventilator circuitry
- **#403 is made specifically for T-Pieces**
- Lifts and supports at the 15 mm site, removing a significant amount of pressure on the distal end of the trach tube – intended to reduce back tracheal wall irritation and initiate cough reflex
- Stabilizes the interface – intended to reduce patient discomfort, movement, and subsequent pulling of the trach, potentially reducing decannulation
- Replaces rubberbands, shoe strings, and other ad hoc securement devices
- **LATEX FREE** - No latex reactions which prevents latex hypersensitivity
- 100% hypoallergenic cotton
- Quick release in case of emergency
- Reimbursement Code A7526

### **DIMENSIONS**

Length: 22" / Width: 1"

### **ORDERING INFORMATION:**

1-800-647-0172 | 610-692-8468 Fax

orders@peppermedical.com | www.PepperMedical.com

# #405

## **VENT-TIE® Adult Trach Tube Holder with Attached Extra Soft Ventilator Anti-Disconnect Strap - 1 Piece**



### Features & Benefits:

- Provides trach tube security while also securing the ventilator circuitry
- Lifts and supports at the 15 mm site, removing a significant amount of pressure on the distal end of the trach tube – intended to reduce back tracheal wall irritation and initiate cough reflex
- Stabilizes the interface – intended to reduce patient discomfort, movement, and subsequent pulling of the trach, potentially reducing decannulation
- Replaces rubberbands, shoe strings, and other ad hoc securement devices
- **LATEX FREE** - No latex reactions which prevents latex hypersensitivity
- 100% hypoallergenic cotton
- Quick release in case of emergency
- Reimbursement Code A7526

### DIMENSIONS

Length: 11.5" / Width: .75"

### ORDERING INFORMATION:

1-800-647-0172 | 610-692-8468 Fax

orders@peppermedical.com | www.PepperMedical.com

# #405W

## **VENT-STRAP® Bariatric Extra Wide/Long Trach Tube Holder with Attached Extra Soft Ventilator Anti-Disconnect Strap - 2 Piece**



### Features & Benefits:

- Provides trach tube security while also securing the ventilator circuitry
- Lifts and supports at the 15 mm site, removing a significant amount of pressure on the distal end of the trach tube – intended to reduce back tracheal wall irritation and initiate cough reflex
- Stabilizes the interface – intended to reduce patient discomfort, movement, and subsequent pulling of the trach, potentially reducing decannulation
- Replaces rubberbands, shoe strings, and other ad hoc securement devices
- **LATEX FREE** - No latex reactions which prevents latex hypersensitivity
- Moisture Repellent Driwick
- Quick release in case of emergency
- Reimbursement Code A7526

### **DIMENSIONS**

Length: 28" / Width: 1.5"

### **ORDERING INFORMATION:**

1-800-647-0172 | 610-692-8468 Fax

orders@peppermedical.com | www.PepperMedical.com

# #406B

## **VENT-TIE II® Adult Extra Wide Trach Tube Holder with Ventilator Anti-Disconnect Device - 2 Piece**



### Features & Benefits:

- Provides trach tube security while also securing the ventilator circuitry
- Lifts and supports at the 15 mm site, removing a significant amount of pressure on the distal end of the trach tube – intended to reduce back tracheal wall irritation and initiate cough reflex
- Stabilizes the interface – intended to reduce patient discomfort, movement, and subsequent pulling of the trach, potentially reducing decannulation
- Replaces rubberbands, shoe strings, and other ad hoc securement devices
- **LATEX FREE** - No latex reactions which prevents latex hypersensitivity
- 100% hypoallergenic cotton
- Quick release in case of emergency
- Reimbursement Code A7526

### **DIMENSIONS**

Length: 22" / Width: 1.5"

### **ORDERING INFORMATION:**

1-800-647-0172 | 610-692-8468 Fax

orders@peppermedical.com | www.PepperMedical.com

# #497

## VENT-TIE® Adult With 1 Piece Trach Tie



### Features & Benefits:

- Provides trach tube security while also securing the ventilator circuitry
- Lifts and supports at the 15 mm site, removing a significant amount of pressure on the distal end of the trach tube – intended to reduce back tracheal wall irritation and initiate cough reflex
- Stabilizes the interface – intended to reduce patient discomfort, movement, and subsequent pulling of the trach, potentially reducing decannulation
- Replaces rubberbands, shoe strings, and other ad hoc securement devices
- **LATEX FREE** - No latex reactions which prevents latex hypersensitivity
- 0.25" Thick blue foam with moisture repellent backing
- Quick release in case of emergency
- Reimbursement Code A7526

### DIMENSIONS

497S: Length: 9.5" / Width: 1.25"

497M: Length: 18.5" / Width: 1.25"

497L: Length: 23.5" / Width: 1.25"

### ORDERING INFORMATION:

1-800-647-0172 | 610-692-8468 Fax

orders@peppermedical.com | www.PepperMedical.com

# #601

## **FOLEY SWIVEL LOC®** **Foley Catheter** **Stabilization Device**



### Features & Benefits:

- Swivel locking mechanism allows for rotation of the Foley Catheter
- Locking system reduces potential of catheter kinking and twisting
- Minimizes dislodgements, traction and positioning
- **LATEX FREE** - No latex reactions which prevents latex hypersensitivity
- Addresses CDC's guidelines for securing Foley Catheters
- Reimbursement Code A4334

### **DIMENSIONS**

Length: 3" / Width: 2"

### **ORDERING INFORMATION:**

1-800-647-0172 | 610-692-8468 Fax

orders@peppermedical.com | www.PepperMedical.com

# #606

## FOLEY-TIE® Foley Catheter Legband



### Features & Benefits:

- One piece Velcro® blue locking tab “Flex Point”
- Soft stretch leg band material
- Quick and easy application around leg and Foley
- No adhesive material
- **LATEX FREE** - No latex reactions which prevents latex hypersensitivity
- Fits up to 24” (most adults and children)
- Reimbursement Code A4334

#### DIMENSIONS

Length: 35” / Width: 2”

#### ORDERING INFORMATION:

1-800-647-0172 | 610-692-8468 Fax

orders@peppermedical.com | www.PepperMedical.com

# #606XL

## FOLEY-TIE® Foley Catheter Bariatric Legband



### Features & Benefits:

- One piece Velcro® blue locking tab “Flex Point”
- Soft stretch leg band material
- Quick and easy application around leg and Foley
- No adhesive material
- **LATEX FREE** - No latex reactions which prevents latex hypersensitivity
- Fits up to 60”
- Reimbursement Code A4334

#### DIMENSIONS

Length: 87.5” / Width: 2”

#### ORDERING INFORMATION:

1-800-647-0172 | 610-692-8468 Fax

orders@peppermedical.com | www.PepperMedical.com

# #5445

## TUBE-TIE® Patch – Multi-Purpose Adhesive Tube Holder With Quick Release



### Applications Include:

- Foley Catheter holder
- Jejunostomy & Gastrostomy tube holders
- Wound drainage tube
- Subclavian line & other ports

### Features & Benefits:

- Securely locks catheter in place
- Non-slip Velcro® Locking Tab can be repositioned
- Secures all types of tubing
- Produces quick and effective removals
- **LATEX FREE** - No latex reactions which prevents latex hypersensitivity
- Individually packaged in boxes of 50 each
- Reimbursement Code A4334

#### DIMENSIONS

Length: 3" / Width: 3"


#### ORDERING INFORMATION:

1-800-647-0172 | 610-692-8468 Fax

orders@peppermedical.com | www.PepperMedical.com

# #650

## Foley Loc® Adhesive Foley Catheter Holder



### Features & Benefits:

- Kink Resistant – non-swivel locking system reduces potential of catheter kinking and twisting
- Adjustable – Velcro® EZ-Tabs allow for release, inspection / cleaning and re-attachment
- Security – minimizes dislodgements, traction and pistoning
- Versatility – secures and stabilizes all styles of Foley catheters (silicone or latex, two or three port catheters)
- Universal Design – one style reduces inventory and save costs
- **LATEX FREE** and hypo-allergenic medical grade adhesive
- Addresses CDC's guidelines for securing Foley Catheters
- Individually packaged in boxes of 50 each
- Reimbursement Code A4333

### DIMENSIONS

Length: 1.75" / Width: 5.5"

### ORDERING INFORMATION:

1-800-647-0172 | 610-692-8468 Fax

orders@peppermedical.com | www.PepperMedical.com

# #5447

## TUBE LOC® Universal Tube Holder



### Features & Benefits:

- Secures any of the following tubes: PICC Lines, Central Venous Catheters, IV Lines, Foley Catheters, Feeding Tubes, Surgical Drains, and Chest Tubes
- Minimizes dislodgements, traction and positioning
- **LATEX FREE** - No latex reactions which prevents latex hypersensitivity
- Addresses CDC's guidelines for tube securement
- Reimbursement Code A4333

#### DIMENSIONS

Length: 2.5" / Width: 1.5"

#### ORDERING INFORMATION:

1-800-647-0172 | 610-692-8468 Fax

orders@peppermedical.com | www.PepperMedical.com

# HME's & Filters

## Heat Moisture Exchangers and Filters

ITEM#	DESCRIPTION
HME-01P	Heat Moisture Exchanger for Trach, Paper
HME-02P	Heat Moisture Exchanger for Trach, Paper with O2 Port
HME-02-Sux	Heat Moisture Exchanger for Trach, Paper with O2 and Suction Port
HME-01F	Heat Moisture Exchanger for Trach, Foam
HME-02F	Heat Moisture Exchanger for Trach, Foam with O2 Port
HME-03P	Heat Moisture Exchanger for Trach, with O2 and Suction Port
HME-04P	Heat Moisture Exchanger for Trach, with O2 and Suction Port
HME-Mini-P	Pediatric/Neo Heat Moisture Exchanger for Trach, Paper
HEPA701	HEPA Filter, In Line
HME-711-Mini	In Line Heat Moisture Exchanger
HME-712	In Line Heat Moisture Exchanger with Filter
HME-712-Mini	In Line Heat Moisture Exchanger, Pediatric
HME-712-Flex	In Line Heat Moisture Exchanger with Flex Tube
BV001	Bacterial Viral Filter
BV002	Bacterial Viral Filter with O2 Port

### ORDERING INFORMATION:

1-800-647-0172 | 610-692-8468 Fax

orders@peppermedical.com | www.PepperMedical.com

# HME-01P

## Heat Moisture Exchanger with Paper Filter



### Description:

Product Code:	<b>HME-01P</b>
Description:	<b>Tracheal HME</b>
Humidity Output:	<b>23.4mg @ VT 500ml</b>
Resistance @ 30 LPM:	<b>40pa</b>
Dead Space:	<b>11ml</b>
Tidal Volume Range:	<b>&gt; 25ml</b>

### Features & Benefits:

- Low profile and lightweight
- High performance paper element
- Warms and humidifies air with minimal resistance to flow
- Convenient oxygen delivery
- Full compatibility with breathing systems

#### ORDERING INFORMATION:

1-800-647-0172 | 610-692-8468 Fax  
orders@peppermedical.com | www.PepperMedical.com

# HME-02P

## Heat Moisture Exchanger with Paper Filter and O2 Port



### Description:

Product Code:	<b>HME-02P</b>
Description:	<b>Tracheal HME</b>
Humidity Output:	<b>23.4mg @ VT 500ml</b>
Resistance @ 30 LPM:	<b>40pa</b>
Dead Space:	<b>11ml</b>
Tidal Volume Range:	<b>&gt; 25ml</b>

### Features & Benefits:

- Low profile and lightweight
- High performance paper element
- Warms and humidifies air with minimal resistance to flow
- **Integrated oxygen port**
- Full compatibility with breathing systems

#### ORDERING INFORMATION:

1-800-647-0172 | 610-692-8468 Fax  
orders@peppermedical.com | www.PepperMedical.com

# HME-02-Sux

## Heat Moisture Exchanger with Paper Filter, O2 Port and Suction Port



### Description:

Product Code:	<b>HME-02-Sux</b>
Dimensions:	<b>H: 38.5mm / W: 28.5mm</b>
Humidity Output:	<b>23.4mg @ VT 500ml</b>
Resistance @ 30 LPM:	<b>40pa</b>
Tidal Volume Range:	<b>&gt; 25ml</b>
Media:	<b>Corrugated paper</b>
Weight:	<b>5g</b>
Connection:	<b>15mm</b>

### Features & Benefits:

- Low profile and lightweight
- High performance paper element
- Warms and humidifies air with minimal resistance to flow
- **Integrated oxygen and suction ports**
- Full compatibility with breathing systems

#### ORDERING INFORMATION:

1-800-647-0172 | 610-692-8468 Fax  
orders@peppermedical.com | www.PepperMedical.com

# HME-01F

## Heat Moisture Exchanger with Foam Filter



### Description:

Product Code:	<b>HME-01F</b>
Description:	<b>Tracheal HME</b>
Humidity Output:	<b>23.4mg @ VT 500ml</b>
Resistance @ 30 LPM:	<b>40pa</b>
Dead Space:	<b>11ml</b>
Tidal Volume Range:	<b>&gt; 25ml</b>
Weight:	<b>5g</b>

### Features & Benefits:

- Low profile and lightweight
- High performance paper element
- Warms and humidifies air with minimal resistance to flow
- Convenient oxygen delivery
- Full compatibility with breathing systems

#### ORDERING INFORMATION:

1-800-647-0172 | 610-692-8468 Fax

orders@peppermedical.com | www.PepperMedical.com

# HME-02F

## Heat Moisture Exchanger with Foam Filter and O2 Port



### Description:

Product Code:	<b>HME-02F</b>
Description:	<b>Tracheal HME</b>
Humidity Output:	<b>23.4mg @ VT 500ml</b>
Resistance @ 30 LPM:	<b>40pa</b>
Dead Space:	<b>11ml</b>
Tidal Volume Range:	<b>&gt; 25ml</b>
Weight:	<b>5g</b>

### Features & Benefits:

- Low profile and lightweight
- High performance paper element
- Warms and humidifies air with minimal resistance to flow
- **Integrated oxygen port**
- Full compatibility with breathing systems


#### ORDERING INFORMATION:

1-800-647-0172 | 610-692-8468 Fax

orders@peppermedical.com | www.PepperMedical.com

# HME-03P

## Heat Moisture Exchanger with Fold Away O2 Port and Suction Port



### Description:

Product Code:	<b>HME-03P</b>
Description:	<b>Tracheal HME</b>
Humidity Output:	<b>23.4mg @ VT 500ml</b>
Resistance @ 30 LPM:	<b>40pa</b>
Dead Space:	<b>11ml</b>
Tidal Volume Range:	<b>&gt; 25ml</b>
Weight:	<b>5g</b>

### Features & Benefits:

- Low profile and lightweight
- High performance paper element
- Warms and humidifies air with minimal resistance to flow
- **Fold away oxygen port plus suction port**
- Full compatibility with breathing systems

#### ORDERING INFORMATION:

1-800-647-0172 | 610-692-8468 Fax

orders@peppermedical.com | www.PepperMedical.com

# HME-04P

## Heat Moisture Exchanger with Side O2 Port and Suction Port



### Description:

Product Code:	<b>HME-04P</b>
Description:	<b>Tracheal HME</b>
Humidity Output:	<b>23.4mg @ VT 500ml</b>
Resistance @ 30 LPM:	<b>40pa</b>
Dead Space:	<b>11ml</b>
Tidal Volume Range:	<b>&gt; 25ml</b>
Weight:	<b>5g</b>

### Features & Benefits:

- Low profile and lightweight
- High performance paper element
- Warms and humidifies air with minimal resistance to flow
- **Side oxygen port plus suction port**
- Full compatibility with breathing systems

#### ORDERING INFORMATION:

1-800-647-0172 | 610-692-8468 Fax

orders@peppermedical.com | www.PepperMedical.com

# HME-Mini-P

## Neonatal and Pediatric Heat Moisture Exchanger



### Description:

Product Code:	<b>HME-Mini-P</b>
Description:	<b>Mini HME</b>
Humidity Output:	<b>32mg H2O @ VT 25ml</b>
Resistance @ 5 LPM:	<b>90pa</b>
Dead Space:	<b>3ml</b>
Tidal Volume Range:	<b>&gt; 10ml</b>
Weight:	<b>3g</b>

### Features & Benefits:

- Low profile and lightweight
- Low Volume Tracheal HME
- Connects to a 15F-15M connection
- **Specifically for Neonatal and Pediatric Patients**

#### ORDERING INFORMATION:

1-800-647-0172 | 610-692-8468 Fax  
orders@peppermedical.com | www.PepperMedical.com

# HEPA701

## HEPA Filter with HME



### Description:

Product Code:	<b>HEPA701</b>
Description:	<b>Pleated HEPA Filter</b>
Humidity Output:	<b>25mg H2O @ VT 500ml</b>
Resistance @ 30 LPM:	<b>140pa</b>
Dead Space:	<b>55ml</b>
Tidal Volume Range:	<b>150-1500ml</b>
Filtration Efficiency:	<b>BFE 99.999% / VFE 99.999%</b>
Weight:	<b>40g</b>
Dimensions:	<b>H: 81.5mm / W: 68.5mm</b>
Materials of Construction:	<b>Filter Media: Hydrophobic Media Housing: Polypropylene</b>
Gas Monitoring Luer Port:	<b>Yes</b>
Connections:	<b>22M-15F / 22F-15M</b>

The Pepper Medical HEPA Filter come with efficiencies of 99.99%, with luer ports, and an integral cap-and-strap as an additional safety feature to remove the risk of a lost cap. Pleated media filters provide higher efficiencies and greater levels of hydrophobicity.

#### ORDERING INFORMATION:

1-800-647-0172 | 610-692-8468 Fax

orders@peppermedical.com | www.PepperMedical.com

# HME711-Mini

## Mini Heat Moisture Exchanger



### Description:

Product Code:	<b>HME711-Mini</b>
Description:	<b>HME</b>
Humidity Output:	<b>10.5mg @ VT 500ml</b>
Resistance to Flow:	<b>40pa @ 30 LPM</b>
Dead Space:	<b>38ml</b>
Tidal Volume Range:	<b>75-1000ml</b>
Filtration:	<b>None</b>

### ORDERING INFORMATION:

1-800-647-0172 | 610-692-8468 Fax  
orders@peppermedical.com | www.PepperMedical.com

# HME-712

## In-Line Heat Moisture Exchanger



### Description:

Product Code:	<b>HME-712</b>
Description:	<b>Paper HME</b>
Humidity Output:	<b>10 @ VT 500ml</b>
Resistance to Flow:	<b>145 @ 30 LPM</b>
Dead Space:	<b>52ml</b>
Tidal Volume Range:	<b>150-1500ml</b>
Materials of Construction:	<b>Paper</b>
Gas Port:	<b>Yes</b>
Connections:	<b>22M-15F / 22F-15M</b>

### ORDERING INFORMATION:

1-800-647-0172 | 610-692-8468 Fax  
orders@peppermedical.com | www.PepperMedical.com

# HME-712-Mini

## In-Line Heat Moisture Exchanger



### Description:

Product Code	<b>#HME-712-Mini</b>
Description	<b>Paper HME, In-line</b>
Humidity Output:	<b>27mg(@Vt 500ml)</b>
Resistance to Flow :	<b>300pa at 60 l/min</b>
Tidal Volume	<b>&lt;600ml</b>
Media:	<b>Paper</b>
Gas Port	<b>No</b>
Dead Space	<b>11ml</b>
Connection:	<b>22M-15F / 22F-15M</b>

### ORDERING INFORMATION:

1-800-647-0172 | 610-692-8468 Fax  
orders@peppermedical.com | www.PepperMedical.com

# HME-712-Flex

## In-Line Heat Moisture with BuiltIn Flex Tube



### Description:

Product Code:	<b>HME-712-Flex</b>
Description:	<b>HME with Flex Tubing</b>
Humidity Output:	<b>4.2 @ VT 500ml</b>
Resistance @ 30 LPM:	<b>20 @ 30 LPM</b>
Tidal Volume Range:	<b>150-1500ml</b>

### ORDERING INFORMATION:

1-800-647-0172 | 610-692-8468 Fax  
orders@peppermedical.com | www.PepperMedical.com

# BV001

## Bacterial Viral Filter



### Description:

Product Code:	<b>BV001</b>
Description:	<b>Bacterial Viral Filter</b>
Weight:	<b>G.W.: 5 / N.W.: 4</b>
Filtration Efficiency:	<b>BFE 99.9999% / VFE 99.9999%</b>
Gas Port:	<b>Yes</b>
Connections:	<b>22M-15F / 22F-15M</b>

### ORDERING INFORMATION:

1-800-647-0172 | 610-692-8468 Fax  
orders@peppermedical.com | www.PepperMedical.com

# BV002

## Bacterial Viral Filter



### Description:

Product Code:	<b>BV002</b>
Description:	<b>Bacterial Viral Filter</b>
Weight:	<b>G.W.: 5 / N.W.: 4</b>
Filtration Efficiency:	<b>BFE 99.99% / VFE 99.99%</b>
Gas Port:	<b>Yes</b>
Connections:	<b>22M-15F / 22F-15M</b>

### ORDERING INFORMATION:

1-800-647-0172 | 610-692-8468 Fax  
orders@peppermedical.com | www.PepperMedical.com

# RESPIRATORY MASKS



Trach Masks



Aerosol Masks



Venturi Masks

Oxygen Masks

## Varieties Available:

- **Aerosol Masks** - 360 degree swivel connector for maximum patient comfort
- **Venturi Masks** - Rotating orifice which is able to set up the flow rates and percentage concentration required
- **Oxygen Masks** - Swivel oxygen mask with CO2 sampling line and filter
- **Trach Masks** - 360 degree swivel connector with 6mm connector

## Features and Benefits:

- 100% medical grade material
- Kink-free, crush-resistant tubing
- Clear, soft vinyl for patient comfort and visual patient assessment
- Material: DEHP free, latex-free

### ORDERING INFORMATION:

1-800-647-0172 | 610-692-8468 Fax  
orders@peppermedical.com | www.PepperMedical.com

# OXYGEN TUBING

---



## Features and Benefits:


- 100% medical grade material
- Kink-free, crush-resistant tubing
- Material: DEHP free, latex-free
- Universal connectors on both sides
- Available in various sizes (7 ft, 25 ft, etc.)
- Available with and without cannulas
- Available for adult, pediatric and infant
- Available with CO2 sampling line

### ORDERING INFORMATION:


1-800-647-0172 | 610-692-8468 Fax  
orders@peppermedical.com | www.PepperMedical.com

# BREATHING CIRCUIT TUBING

---



**Corrugated  
Taper Tubing**



**Expandable  
Taper Tubing**

## Benefits:

- Tubing for breathing circuit systems
- Fits standard 15 mm and 22 mm connections
- Available in two configurations:

---

**Circuit Extension | Expandable, 6"**

---

**Circuit Extension | Corrugated, 6"**

### ORDERING INFORMATION:

1-800-647-0172 | 610-692-8468 Fax

orders@peppermedical.com | www.PepperMedical.com

# #654

## PULSE OXIMETER FINGER WRAP



### Benefits:

- Pulse oximeter probe wrap that helps secure various styles of optical transducers and probe sensors
- Prevents the misalignment of optical components while blocking out extra light
- Fits a neonatal foot, pediatric toes, or pediatric and adult fingers
- Packaged 12 wraps per package
- Velcro® tab eliminates need for adhesive tape
- Skin friendly foam laminate made from moisture repellent material
- Over-the-counter intended for non-medical, recreational use
- Adjustable, one size fits all
- Secures Pulse Ox sensor to allow for steady reading
- Measures 5" long x 1" wide and can be cut to fit any length

### ORDERING INFORMATION:

1-800-647-0172 | 610-692-8468 Fax

orders@peppermedical.com | www.PepperMedical.com

# TRACHEOSTOMY SPEAKING VALVES

The SPIRO Speaking Valve is a lightweight, low profile, one-way valve designed to be attached to the 15mm connector of the tracheostomy tube. Optional Oxygen port and particle filters are also available.

## #710

### SPEAKING VALVE\*\*

#### Features & Benefits:

- Speaking Valve with oxygen port for spontaneous breathing tracheostomy patients.
- Intended to allow the patient to speak without the need for finger occlusion. Also provides back pressure for other common body functions. Exchangeable filter to prevent course particles entering the respiratory system.
- This transparent device is recommended for use in the home care environment.

*\*\*Also available without O2 Port*



## #711

### SPEAKING VALVE

#### Features & Benefits:

- Speaking Valve for spontaneously breathing tracheostomy patients.



## #740

### REPLACEMENT FILTER (OPTIONAL)

#### Features & Benefits:

- Exchangeable filter which is compatible 710, 710R.
- The filter prevents course particles from entering respiratory system.





**ORDERING INFORMATION:** 1-800-647-0172 | 610-692-8468 Fax  
orders@peppermedical.com | www.PepperMedical.com

Manufactured by Fogless International AB for Pepper Medical


# UNI-BLADE

*One Blade Fits All 4 Sizes*


- ▶ *Eliminates the need for multiple laryngoscope blades*
- ▶ *Saves space, time and money*




Size 1



Size 2



Size 3



Size 4

## Features and Benefits:

- **Patented One-Blade-Fits-All Design:** Uni-Blades unique multi-positional attachment system allows for functional adjustment of the blade length.
- **Eliminate the Need for Multiple Laryngoscope Blades:** This system ensures a precise, one-blade-fits-all solution.
- **One Adjustable Blade with 4 Size Positions:** Uni-Blades adjustable length ensures the right fit for each patient, streamlining the process in anesthesia, emergency and clinical settings.
- **Space Saving Design:** By combining 4 blade sizes into one, the Uni-Blade reduces clutter on crash carts, field kits, ambulances, helicopters, and other transport vehicles. Ideal for military medical units where space and weight are critical. Ideal solution for any inventory management system. Uni-Blade reduces the number of inventory items stocked reducing cost, time and maximizing inventory efficiencies.
- **High-Strength Curved Material with LED Light Tube:** Combines durability and superior visibility to improve patient outcomes during intubation.
- **Disposable Design:** Ensures hygiene and reduces the risk of cross-contamination in high-risk environments.

[Click Here to Watch the Uni-Blade Video](#)


### ORDERING INFORMATION:

1-800-647-0172 | 610-692-8468 Fax

orders@peppermedical.com | www.PepperMedical.com

## #701VCS #701VCS-Pediatric

# VENTILATOR CIRCUIT STABILIZER



## Features & Benefits:

- The Ventilator Circuit Stabilizer secures the vent tubing to the patients' chest or abdomen, providing for stabilization of the tubing that won't allow it to become displaced during patient care or during routine patient movement.
- The VCS is a belt made of a soft cotton laminate that fits comfortably around a patient's waist.
- It comes with a strap that is permanently fastened on one end to the belt while the other end of the strap is fitted with a small piece of VELCRO® brand hook.
- The vent circuit is secured to the belt by laying the circuit onto the belt and then placing the strap over the circuit and attaching the hook end of the strap to the belt.
- Also available with an attached nasogastric tube securement strap.
- Patent pending


### ORDERING INFORMATION:

1-800-647-0172 | 610-692-8468 Fax

orders@peppermedical.com | www.PepperMedical.com

## #701VCS/NG

# VENTILATOR CIRCUIT AND NG TUBING STABILIZER (Pediatric)



### Features & Benefits:

- The VCS/NG is a belt made of a soft cotton laminate that fits comfortably around a patient's waist.
- It comes with two straps that are permanently fastened on one end to the belt while the other ends of the strap are fitted with a small piece of VELCRO® brand hook.
- The Vent Circuit Tubing is secured to the belt by laying the circuit onto the belt and then placing the strap over the circuit and attaching the hook end of the strap to the belt.
- The NasoGastric Tube is secured to the belt by laying the circuit onto the belt and then placing the strap over the circuit and attaching the hook end of the strap to the belt.
- Simple to use, efficient, and comfortable
- Patent pending

#### ORDERING INFORMATION:

1-800-647-0172 | 610-692-8468 Fax

orders@peppermedical.com | www.PepperMedical.com

# 3 AND 4 PANEL ABDOMINAL BINDERS



## Features & Benefits:

- LATEX FREE
- Comfortable, stretchable supports
- Aids in redevelopment of abdominal muscles
- Evenly distributed compression encourages deep breathing
- Reduces complications such as pneumonia, phlebitis, and atelectasis
- Easy to use fasteners
- Reimbursement Code A4462

ITEM #	DESCRIPTION	DIMENSIONS	MATERIAL	CODE
#AB330	3 Panel Ab Binder, Small	Length: 30" - 45" Width: 9"	Elastic and hook/loop	A4462
#AB345	3 Panel Ab Binder, Medium	Length: 45" - 60" Width: 9"	Elastic and hook/loop	A4462
#AB360	3 Panel Ab Binder, Large	Length: 60" - 75" Width: 9"	Elastic and hook/loop	A4462
#AB430	4 Panel Ab Binder, Small	Length: 30" - 45" Width: 12"	Elastic and hook/loop	A4462
#AB445	4 Panel Ab Binder, Medium	Length: 45" - 60" Width: 12"	Elastic and hook/loop	A4462
#AB460	4 Pabel Ab Binder, Large	Length: 60" - 75" Width: 12"	Elastic and hook/loop	A4462

### ORDERING INFORMATION:

1-800-647-0172 | 610-692-8468 Fax

orders@peppermedical.com | www.PepperMedical.com

# ABDOMINAL BELTS



## Features & Benefits:

- LATEX FREE
- Comfortable, stretchable supports
- Aid in redevelopment of abdominal muscles
- Evenly distributed compression encourages deep breathing
- Reduces complications such as pneumonia, phlebitis, and atelectasis
- Easy to use fasteners
- Reimbursement Code A4462

ITEM #	DESCRIPTION	DIMENSIONS	MATERIAL	CODE
#RB2F	Rib Belt, 2 Panel, Female	Length: 28" - 50" 1 Lb.	Elastic and hook/loop	A4462
#RB2M	Rib Belt, 2 Panel, Male	Length: 45" - 60" 1 Lb.	Elastic and hook/loop	A4462

### ORDERING INFORMATION:

1-800-647-0172 | 610-692-8468 Fax

orders@peppermedical.com | www.PepperMedical.com