

TRACH-TIE®

Tracheostomy Tube Neckband Catalog



PEPPER MEDICAL

1-800-647-0172 | PepperMedical.com

Contents

Page 2 | #**300**

Pages 3 - 4 | #**301 & #301P**

Pages 5 - 7 | #**501, #501P & #501W**

Pages 8 - 9 | #**502 & #502P**

Page 10 | #**503 & 504**

Page 11 | #**540, #541 & #542**

Page 12 | #**597S, #597M & #597L**

Page 13 | #**1001**

Page 14 | #**301PRV & #501PRV**

#300

TRACH-TIE® **Adult Tracheostomy Tube** **Neckband - 1 Piece**



Features & Benefits:

- **LATEX FREE**
- Soft, moisture repellent, fabric ensures patient safety and comfort
- Soft, plush backing on hook tabs provides additional comfort
- Reduces skin breakdown
- Easy to adjust for proper fit up to a 21" neck
- "Closed Seams" on edge of neckband reduces sponge effect of foam laminate and keeps moisture out of trach tie
- Thicker foam padding
- Reimbursement Code A7526

DIMENSIONS

Length: 17" / Width: 1"

ORDERING INFORMATION:

1-800-647-0172 | 610-692-8468 Fax

orders@peppermedical.com | www.PepperMedical.com

#301

TRACH-TIE® **Adult Tracheostomy Tube** **Neckband - 1 Piece**



Features & Benefits:

- **LATEX FREE**
- Soft, 100% hypoallergenic, cotton fabric ensures patient safety and comfort
- Reduces skin breakdown
- Easy to adjust for proper fit up to a 21" neck
- "Closed Seams" on edge of neckband reduces sponge effect of foam laminate and keeps moisture out of trach tie
- Thicker foam padding
- Reimbursement Code A7526

DIMENSIONS

Length: 17" / Width: 1"

ORDERING INFORMATION:

1-800-647-0172 | 610-692-8468 Fax

orders@peppermedical.com | www.PepperMedical.com

#301P

PEDI-TIE®

Pediatric Tracheostomy Tube Neckband - 1 Piece



Features & Benefits:

- **LATEX FREE**
- Soft, 100% hypoallergenic, cotton fabric ensures patient safety and comfort
- Reduces skin breakdown
- **#301P is made specifically for the Pediatric Patient**
- “Closed Seams” on edge of neckband reduces sponge effect of foam laminate and keeps moisture out of trach tie
- Thicker foam padding
- Reimbursement Code A7526

DIMENSIONS

Length: 11.5” / Width: .75”

ORDERING INFORMATION:

1-800-647-0172 | 610-692-8468 Fax

orders@peppermedical.com | www.PepperMedical.com

#501

TRACH-TIE® **Adult Tracheostomy Tube** **Neckband - 2 Piece**



Features & Benefits:

- **LATEX FREE**
- Soft, 100% hypoallergenic, cotton fabric ensures patient safety and comfort
- Reduces skin breakdown
- Two piece adjustable neckband for universal sizing
- “Closed Seams” on edge of neckband reduces sponge effect of foam laminate and keeps moisture out of trach tie
- “Flex-Piece” for cough reflex
- Thicker foam padding
- Reimbursement Code A7526

DIMENSIONS

Length: 22.5” / Width: 1”

ORDERING INFORMATION:

1-800-647-0172 | 610-692-8468 Fax
orders@peppermedical.com | www.PepperMedical.com

#501P

PEDI-TIE II® **Pediatric Tracheostomy Tube** **Neckband - 2 Piece**



Features & Benefits:

- **LATEX FREE**
- Soft, 100% hypoallergenic, cotton fabric ensures patient safety and comfort
- Reduces skin breakdown
- **#501P is made specifically for the Pediatric Patient**
- Two piece adjustable neckband for universal sizing
- “Closed Seams” on edge of neckband reduces sponge effect of foam laminate and keeps moisture out of trach tie
- Thicker foam padding
- Reimbursement Code A7526

DIMENSIONS

Length: 14” / Width: .75”

ORDERING INFORMATION:

1-800-647-0172 | 610-692-8468 Fax
orders@peppermedical.com | www.PepperMedical.com

#501W

TRACH-TIE® Bariatric Tracheostomy Tube Neckband - 2 Piece



Features & Benefits:

- **LATEX FREE**
- 1.5" WIDE design prevents roping effect
- Two piece adjustable neckband allows for easy adjustment
- Soft, plush backing on hook tabs provides additional comfort
- Made of DriWick®, a proprietary material that pulls moisture away for the patients skin, reducing skin breakdown and irritation
- Extra long for bariatric patients
- Reimbursement Code A7526

DIMENSIONS

Length: 28" / Width: 1.5"

ORDERING INFORMATION:

1-800-647-0172 | 610-692-8468 Fax

orders@peppermedical.com | www.PepperMedical.com

#502

TRACH-TIE® **Adult Tracheostomy Tube** **Neckband - 2 Piece**



Features & Benefits:

- **LATEX FREE**
- Soft, 100% hypoallergenic, cotton fabric ensures patient safety and comfort
- 2.25" Velcro® hook tabs, tapered to prevent skin breakdown
- Two piece adjustable neckband for universal sizing
- 0.25" thick foam with moisture repellent backing
- Reimbursement Code A7526

DIMENSIONS

Length: 18.5" / Width: 1"

ORDERING INFORMATION:

1-800-647-0172 | 610-692-8468 Fax

orders@peppermedical.com | www.PepperMedical.com

#502P

PEDI-TIE® **Pediatric Tracheostomy Tube** **Neckband - 2 Piece**



Features & Benefits:

- **LATEX FREE**
- Soft, 100% hypoallergenic, cotton fabric ensures patient safety and comfort
- 2.25" Velcro® hook tabs, tapered to prevent skin breakdown
- Two piece adjustable neckband for universal sizing
- 0.25" thick foam with moisture repellent backing
- Reimbursement Code A7526

DIMENSIONS

Length: 11" / Width: .75"


ORDERING INFORMATION:

1-800-647-0172 | 610-692-8468 Fax

orders@peppermedical.com | www.PepperMedical.com

#503

TRACH-TIE® Pediatric Tracheostomy Tube Neckband - 2 Piece



Features & Benefits:

- **LATEX FREE**
- 0.25" thick white foam with 100% cotton laminate
- Reduces skin breakdown
- **#503 is made specifically for the Pediatric Patient**
- Two piece adjustable neckband for universal sizing
- 2 .875" Velcro® hook tabs, stitched
- Reimbursement Code A7526

DIMENSIONS: Length: 12" / Width: 1"

#504

TRACH-TIE® Adult Tracheostomy Tube Neckband - 2 Piece



Features & Benefits:

- **LATEX FREE**
- 0.25" thick white foam with 100% cotton laminate
- Reduces skin breakdown
- Two piece adjustable neckband for universal sizing
- 2 .875" Velcro® hook tabs, stitched
- Reimbursement Code A7526

DIMENSIONS: Length: 18.5" / Width: 1"

ORDERING INFORMATION: 1-800-647-0172 | 610-692-8468 Fax
orders@peppermedical.com | www.PepperMedical.com


#540

TRACH-TIE® Adult Tracheostomy Tube Neckband - 2 Piece

Features & Benefits:

- **LATEX FREE**
- 0.25" thick blue foam with moisture repellent backing
- 2 .875" Velcro® hook tabs, stitched
- Reimbursement Code A7526

DIMENSIONS: Length: 18.5" (including 2.5" of elastic) / Width: 1"



#541

TRACH-TIE® PEDIATRIC Tracheostomy Tube Neckband - 2 Piece

Features & Benefits:

- **LATEX FREE**
- 0.1875" thick blue foam with moisture repellent backing
- 2 .875" Velcro® hook tabs, welded
- Reimbursement Code A7526

DIMENSIONS: Length: 17.5" / Width: 1"



#542

TRACH-TIE® PEDIATRIC / NEONATAL Tracheostomy Tube Neckband - 2 Piece

Features & Benefits:

- **LATEX FREE**
- 0.1875" thick blue foam with moisture repellent backing
- 2 .875" Velcro® hook tabs, welded
- Reimbursement Code A7526

DIMENSIONS: Length: 10" / Width: 0.75"



ORDERING INFORMATION: 1-800-647-0172 | 610-692-8468 Fax
orders@peppermedical.com | www.PepperMedical.com


#597S

TRACH-TIE® Small Tracheostomy Tube Neckband - 1 Piece

Features & Benefits:

- **LATEX FREE**
- 0.25" thick blue foam with moisture repellent backing
- 2 .875" x 0.1875 Velcro® hook tabs
- Reimbursement Code A7526

DIMENSIONS: Length: 9.5" / Width: 1.25"




#597M

TRACH-TIE® Medium Tracheostomy Tube Neckband - 1 Piece

Features & Benefits:

- **LATEX FREE**
- 0.25" thick blue foam with moisture repellent backing
- 2 .875" x 0.25 Velcro® hook tabs
- Reimbursement Code A7526

DIMENSIONS: Length: 18.5" / Width: 1.25"




#597L

TRACH-TIE® Large Tracheostomy Tube Neckband - 1 Piece

Features & Benefits:

- **LATEX FREE**
- 0.25" thick blue foam with moisture repellent backing
- 2 .875" x 0.25 Velcro® hook tabs
- Reimbursement Code A7526


DIMENSIONS: Length: 23.5" / Width: 1.25"



ORDERING INFORMATION: 1-800-647-0172 | 610-692-8468 Fax
orders@peppermedical.com | www.PepperMedical.com

#1001

TRACH-TIE® **Adult Tracheostomy Tube** **Neckband - 2 Piece**



Features & Benefits:

- **LATEX FREE**
- 2.875" Velcro® hook tabs, stitched
- Two piece adjustable neckband for universal sizing
- 0.25" thick blue foam with moisture repellent backing
- Reimbursement Code A7526

DIMENSIONS

Length: 18.5" (including 2.5" of elastic) / Width: 1"

ORDERING INFORMATION:

1-800-647-0172 | 610-692-8468 Fax

orders@peppermedical.com | www.PepperMedical.com

#301PRT

RAINBOW TIES® Pediatric/Neonatal Tracheostomy Tube Neckband - 1 Piece



Features & Benefits:

- 100% hypoallergenic, cotton fabric ensures patient safety and comfort
- **Packaged 10 per bag with 3 different color schemes: Variety, Pastel, & Bold**
- **LATEX FREE**
- Made of Orthowick® which pulls moisture away from skin aiding in the reduction of skin breakdown and infection
- Soft Velcro® tabs and plush backing
- Reimbursement Code A7526

DIMENSIONS: Length: 11.5" / Width: 0.75"

#501PRT

RAINBOW TIES® Pediatric/Neonatal Tracheostomy Tube Neckband - 2 Piece



Features & Benefits:

- 100% hypoallergenic, cotton fabric ensures patient safety and comfort
- **Packaged 10 per bag with 3 different color schemes: Variety, Pastel, & Bold**
- **LATEX FREE**
- Made of Orthowick® which pulls moisture away from skin aiding in the reduction of skin breakdown and infection
- Soft Velcro® tabs and plush backing
- Reimbursement Code A7526

DIMENSIONS: Length: 14" / Width: 0.75"

ORDERING INFORMATION: 1-800-647-0172 | 610-692-8468 Fax
orders@peppermedical.com | www.PepperMedical.com