

TRACHEOSTOMY SPEAKING VALVES

The SPIRO Speaking Valve is a lightweight, low profile, one-way valve designed to be attached to the 15mm connector of the tracheostomy tube. Optional Oxygen port and particle filters are also available.

#710 – SPEAKING VALVE*

Speaking Valve with oxygen port for spontaneous breathing tracheostomy patients.

Intended to allow the patient to speak without the need for finger occlusion. Also provides back pressure for other common body functions. Exchangeable filter to prevent coarse particles entering the respiratory system.

This transparent device is recommended for use in the home care environment.



**Also available without O2 Port*

#711 – SPEAKING VALVE

Speaking Valve for spontaneously breathing tracheostomy patients.



#740 – REPLACEMENT FILTER (OPTIONAL)

Exchangeable filter which is compatible 710, 710R. The filter prevents coarse particles from entering respiratory system.

